

FOR PROFESSIONAL USE ONLY



AMINAL Cleaner A

**CARTRIDGE
CONTAINING POWDER
OF ANHYDROUS
SODIUM CARBONATE
– Na_2CO_3**

FOR PROFESSIONAL USE ONLY



AMINAL Cleaner C

**CARTRIDGE
CONTAINING POWDER
OF ANHYDROUS
CITRIC ACID
– $\text{C}_6\text{H}_8\text{O}_7$**

HAEMODIALYSIS

AMINAL Cleaner A

**CARTRIDGE CONTAINING
POWDER OF ANHYDROUS
SODIUM CARBONATE**

– Na₂CO₃

AMINAL Cleaner C

**CARTRIDGE CONTAINING
POWDER OF ANHYDROUS
CITRIC ACID**

– C₆H₈O₇

TECHNICAL DATA - AMINAL Cleaners

Description	AMINAL Cleaner A Cartridge containing powder of anhydrous sodium carbonate – Na ₂ CO ₃	AMINAL Cleaner C Cartridge containing powder of anhydrous citric acid – C ₆ H ₈ O ₇
Composition	Sodium Carbonate (Na ₂ CO ₃)	Citric Acid (C ₆ H ₈ O ₇)
Intended use	On-line production of sodium carbonate solution intended for the disinfection and cleaning of dialysis machines, as well as for the removal of organic deposits from the fluid path of the dialysis machines.	On-line production of citric acid solution intended for the disinfection and cleaning of dialysis machines, as well as for the removal of calcium and magnesium salts from the fluid path of the dialysis machines.
Packaging	13 g	32 g
Compatible with	Gambro dialysis machines	
Shelf life	2 years	
Storage conditions	5°C - 30°C	

The product is in compliance with the requirements of European Pharmacopoeia and USP.



**ALKALOID
SKOPJE**

Health above all

www.alkaloid.com.mk

CE 2803

Version 02 from 05.2023.

PHARMACEUTICAL CHEMICAL COSMETIC INDUSTRY
Blvd. Aleksandar Makedonski 12, Skopje Republic of North Macedonia
Tel.: +389 2 31 04 041; Fax.: +389 2 31 04 041